

Membrane Gas Separation



Membrane Gas Separation				
Input line	4	1600 cc/min – Up to 150 bar		
Mixer	1 set	4 input line		
Vacuum pump	1 set	discharging the lines		
MFC	4	1600 cc/min (Brooks)		
MFM	4	For check the permeate and reiterate flow (Brooks)		
oven	2 set	Max temperature 180 °C	Oven 1	1 or 2 set of HC or LC cell
			Oven 2	1 set of M series module M (10-20-30)
Working/testing pressure	150/250 bar			
Control system	HMI 14" touch screen – PLC Automation			